


OUR UPCOMING EXHIBITIONS:



Perlove

VET1120 Series

Vet Mobile Digital C-arm System

VET1120 Series

Vet Mobile Digital C-arm System

Newly integrated design of the whole machine
Compact and beautiful structure, flexible movement

Digital Workstation with Compact Design

Professional software special for vet use, Smart and compact design with integrated computer workstation, save space and convenient to operate



Hand Controller

Humanized miniature, mobile operator panel design, free control from the host, practical & convenient

Image Intensifier

Toshiba brand image intensifier with reliable quality together with dense grid, enhancing superior image acquisition



Eyes of Doctors



Smart design with flexible movement, provides wider operating space

X-ray tube

Combined X-ray tube assembly with high output power, no cable design, beautiful structure, easy to maintain

Technical Parameters

Items	VET1120 Series		
High frequency inverter power supply	Power:3.5kW	Main inverter frequency: 40 kHz	Power:5.0kW Main inverter frequency: 40 kHz
Dual focus	Small focus:0.3;Large focus:1.5		Small focus:0.3;Large focus:1.5
Automatic fluoroscopy	tube voltage 40 kV~110kV	adjust automatically	tube voltage 40 kV~120kV adjust automatically
	tube current: 0.3mA~4mA	adjust automatically	tube current: 0.3mA~4mA adjust automatically
Manual fluoroscopy	tube voltage 40 kV~110kV	continuous	tube voltage 40 kV~120kV continuous
	tube current: 0.3mA~4mA	continuous	tube current: 0.3mA~4mA continuous
Pulse fluoroscopy	tube voltage 40 kV~110kV	continuous	tube voltage 40 kV~120kV continuous
	tube current: 0.3mA~8mA	continuous	tube current: 0.3mA~30mA continuous
Photography tube voltage, tube current mAs	40kV~110kV	20~63mA 1~125 mAs	40~120kV 25~100mA 1~180mAs